

How to Use the Shower Bag

1. Assemble the shower bag.

- Clip the shoulder strap to the shower bag using the rings.



- Adjust the shoulder strap so it fits snug without pulling on or moving the driveline.
- To attach the clip style belt – slide the belt through the loop on the side of the bag that will be against the body. It can be worn on the right or left depending on strap positioning.



2. Putting on the Shower Bag

- Unclip the top of the shower bag to open.
- Unzip and over the cover.



- Place batteries into the shower bag.

If you have any questions, please call the VAD Office at 612-863-5638 and ask to speak to a VAD Coordinator.

- Slide the system controller into the pocket on the inside cover with the display facing up.



- Close the lid and zip the cover. Make sure the power cables remain inside the bag and the driveline exits through the protective red tab. DO NOT get the driveline caught in the zipper.



- Close the top lid with the snap.



- Hang the shower bag over your head and shoulder so the bag hangs at your side. Clip the belt around your waist to prevent any movement of the bag.

3. To Remove the Shower Bag

- Remove belt and carefully remove shower bag should strap. Place bag on a stable surface. Towel dry the bag before opening.
- Hand bag and allow to air dry.

If you have any questions, please call the VAD Office at 612-863-5638 and ask to speak to a VAD Coordinator.