

How To Clean Your Entere! (ENFit®) Tube

It is important to clean your ENFit tube after each home tube feeding. Follow these steps:

Step 1

- Gather your supplies:
 - glass or cup with clean water
 - syringe
 - gauze pad
 - a clean toothbrush
 - ENFit feeding tube with moat and open cap.
- Wash your hands.
- Rinse the toothbrush with tap water.

Step 2



Fill the syringe with tap water.

Step 3



Plug the center hole of the feeding tube with the toothbrush bristles. Flush the moat with clean tap water.

Step 4



Rotate the toothbrush in the bottom of the moat.

Step 5



Rinse the feeding tube cap with clean tap water.

Step 6

Put the toothbrush bristles into the feeding tube cap. Rotate the brush in the cap.

Step 7



Wipe the feeding tube port and the cap with gauze.

Step 8

Clean all supplies. Let air dry.