

Metered Dose Inhaler With Chamber

A metered dose inhaler with chamber has advantages. It takes away timing issues by letting you breathe at your own pace and


1. Shake the inhaler canister quickly for 5 seconds.



5. Grasp the inhaler in a good position to firmly spray it.



2. Remove the caps from the inhaler and the chamber.



6. Place the chamber 7. Firmly press the canister once.



if it's right for you.

3. Look to make sure

there are not any

foreign objects in


the inhaler or the

chamber.

4. Insert the inhaler mouthpiece into the inhaler adapter of the chamber.

it gets the medicine to your lungs better.

Together, you and your doctor will decide



8. Inhale a full breath slowly through your mouth while keeping your lips sealed around the mouthpiece. A whistling sound means you are breathing too fast.



- **12.** Remove the inhaler and holding chamber. Replace the caps. Store vour inhaler and chamber in a clean, dry location.

(over)

- - your lips. Gently exhale as fully as is comfortable.

mouthpiece in

teeth and close

your mouth between your

- 9. Hold your breath for 10 seconds if you can.
- **10.** Remove the chamber from your mouth and exhale normally.



11. Wait 30 to 60 seconds. Shake the inhaler again and repeat steps 5 through 10 for the number of puffs ordered by your doctor. Note: If you use albuterol. wait 1 to 2 minutes.

Metered Dose Inhaler With Spacers

A metered dose inhaler with spacers has advantages. It is easier to carry and it costs less than an MDI with a chamber.

- **1.** Shake the inhaler canister quickly for 5 seconds.
- **2.** Remove the cap from the inhaler.





5. Grasp the inhaler in a good position to firmly spray it.



6. Place the spacer mouthpiece in your mouth between your teeth and close your lips. Gently exhale as fully as is comfortable.


- Together, you and your doctor will decide if it's right for you.
 - 3. Look to make sure there are not any foreign objects in the inhaler or the spacer.

7. Firmly press the


once.

inhaler canister

4. Insert the inhaler mouthpiece into one end of the spacer.



8. Inhale a full breath slowly through your mouth.





- 9. Hold your breath for 10 seconds if you can.
- **10.** Remove the spacer from your mouth and exhale normally.



- 11. Wait 30 to 60 seconds. Shake the inhaler again and repeat steps 5 through 10 for the number of puffs ordered by your doctor.
 Note: If you use albuterol, wait 1 to 2 minutes.
- 12. Remove the inhaler and spacer. Replace the cap. Store your inhaler and spacer in a clean, dry location.